


Understanding the Chart View Tab

Add or subtract sub-tab charts, refresh the chart, or print the chart

This allows you to select any study from the current patient's history

The current patient and chart date are always on display



Organize your x-rays in different charts, like FMX or Bitewing

Chart View Tab Notes:

1. Single click the images to select them, then double click the last image to send all selected x-rays to the Diagnosis View Tab
2. Double click one image to send it, and the history of that patient's chart position, to the Diagnosis View Tab
 - a. If that position has no history (if it is the first x-ray in that chart position for that patient) then the whole current series goes to the Diagnosis View Tab.
3. Check marks indicate that the image is selected. Single click it again to de-select it.
4. You can select images in the film strip, as well.