


AMRAD MEDICAL



DFMT | FMT


Premier Radiographic Systems

www.AmRadMedical.com

Features & Benefits



DFMT




Fast and Accurate Patient Positioning

Digital Floor Mounted Tubestand

Easy to read LCD shows SID continuous to table-top and table receptor with discrete SID of 40" and 72" to the wallstand

Cross-table views are easily accomplished using the rotational and transverse positioning capabilities

10 ft. floor track allows for off-table and weight bearing exams and is customizable down to 8 ft.



Compatible with most DR Solutions



Easy Patient Transfer



Convenient Lock-Release Pedals



Standing Knees



Ergonomic Design

Elevating Four-way Float Top Table

Flat, four-way, floating top permits easy transfer and positioning of patients

Designed for all types of applications, available in 400lb and 650lb patient capacities

Four-pedal design allows for easy accessibility to table lock release

Collision-protection and fail-safe locks ensure safe operation

Four-way lock release and elevation handle (DFMT)

Options include:

35.5" wide table top

Four-way lock release handle

Non-elevating version available in 400 lb and 650 lb patient capacities

Rotating Cassette Tray

Wallstand

Center mounting or side mounting of Bucky makes this wallstand flexible for any patient application or installation limitation

Ergonomically designed hand control for ease of positioning wall receptor

Side-mounted handgrips stabilize patients during PA exams

Overhead Grip (DFMT)

Remote Collimator Light

Options include:

Dual-sided Lock Release

Rotating Cassette Tray



35.5" Wide Table Top



Economical Alternative FMT



Table Lock Release Handle with Elevation (FMT, optional)



Overhead Grip for J1000 Wallstand (FMT, optional)

System Options (L to R, above)

The AmRad™ Medical DFMT and FMT Systems can be configured to meet the needs of any department. Options include, but are not limited to:

An extra-wide table top for the S222 elevating table to help patients feel more secure

Four-way lock release and elevation handle for easy patient positioning can be located anywhere on the rails of the S223 non-elevating or S222 elevating tables (DFMT, standard; FMT, optional)

The FMT is an economical alternative to the DFMT with the same versatility and functionality as the DFMT

Patient Overhead Grip for use with the J1000 wallstand helps steady patients for upright exams (DFMT, standard; FMT, optional)



7555 N. Caldwell Avenue
Niles, Illinois 60714

Phone: 773-353-4000
Fax: 773-588-3424
Toll Free: 1-800-729-9729



www.AmRadMedical.com



SAI GLOBAL
ISO 13485
Quality