

DELWORKS

WORKSTATION

DELWORKS is a powerful image acquisition and processing software featuring a user-friendly interface that provides sophisticated medical imaging. Designed with complex system automation, DELWORKS aims to simplify the examination process by making the difficult decisions for you.

Its interface features an intuitive, single console workflow designed to simplify every step of the examination procedure. This provides technologists with fast image preview times, advanced image enhancement capabilities and generator control options - all from a large format LED touchscreen display.

Key Highlights

Power Intel® Core™ Desktop

-3.2GHz i5-6500 processor

Full HD Resolution

-23.0" screen display

High Durability

-5H hardened non-scratch glass

Optimal Medical Image Viewing

LED-backlit IPS Panel

-Allowing wide viewing angles that minimize color shift and contrast changes



Technical Specifications


* All specifications subject to change without notice

DELWORKS Workstation Monitor (Eizo MS236WT)

Cabinet color:	Black, Gray
Display:	23.0" Full-HD resolution (1920 x 1080)
Color support:	16.77 million from a palette of 1.06 billion colors
Viewing angles:	178°, 178°
Display size:	20" (50.91 cm) x 11.3" (28.7 cm)
Pixel pitch:	0.270 x 0.270 mm
Brightness:	260 cd/m ²
Contrast ratio:	1000:1
Response time:	16 ms (on/off)
Weight:	14.55 lbs. (6.6 kg) with stand
Dimensions:	21.9" x 14.4" x 8.1"

DELWORKS Workstation Desktop (Dell OptiPlex 7050)

Processor:	3.2GHz Intel® Core™ i5-6500
Memory:	16GB RAM
Hard disk drive:	1TB 7200RPM
Operating system:	Windows 10 Pro (64 bit) English



www.delmedical.com

New York
28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Illinois
241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

DEL MEDICAL

#9000-DW-Workstation 07.19