

# IMAGEWORKS EVA SELECT DIGITAL SENSOR SOLUTION



## The *ImageWorks EVA Select Digital Sensor System*

*has been completely redesigned for quality and performance*

- Market-leading resolution (20+ line pairs per mm)
- Rounded corners for patient comfort
- Razor thin: 5.4 mm thickness
- Integrates with major software platforms
- USB interface works with your existing system




### *Gorgeous Images Plus a Great Price Equals Return on Investment*



*Universal Direct USB connection easily interfaces with your equipment*

### *Need Imaging Software? We Have You Covered*



***EvaSoft Imaging software integrates all of your images into one system.***

***Easy and intuitive to use. No training needed. No mistakes.***

***Has capability to work with all major practice management platforms.***

***Powerful diagnostic capabilities.***

***One-click insurance claims allows you to work faster.***



# IMAGEWORKS EVA SELECT DIGITAL SENSOR SOLUTION

## Technical Specifications:

	Size 2	Size 1
Gray shades:	4096 (12 bits)	4096 (12 bits)
Imager Resolution:	20+ lp/mm	20+ lp/mm
Optical Resolution:	26 lp/mm	26 lp/mm
Pixel Size:	19 microns x 19 microns	19 microns x 19 microns
Active Area:	36 mm (L) x 25.99 mm (W)	30.0 mm (L) x 20.0 mm (W)
Capsule Dimensions:	42.8 mm (L) x 30.5 mm (W) x 5.4 mm (Th)	36.7 mm (L) x 24.4 mm (W) x 5.4 mm (Th)
Image Size:	4.94 MB	3.16 MB



**ImageWorks**

### Who we are:

We know that most dentists are also business owners who focus on getting the best value for their investment. Based in Elmsford, NY for over 35 years, ImageWorks has been manufacturing dental imaging solutions that provide great return on investment. This value stems not only from the quality of our products, but also the quality of our people.

## The ImageWorks Full Dental Imaging End to End Solution



Imaging Software




Digital Intraoral Sensors



Intraoral X-Rays



Panoramic and  
Cephalometric X-Rays



CBCT Systems



**ImageWorks**

Call 914-592-6100 or email [Custserv@imageworkscorporation.com](mailto:Custserv@imageworkscorporation.com)

Visit us at [www.imageworkscorporation.com](http://www.imageworkscorporation.com)